




Code No.KAL-KR103

For research use only

Anti Mouse CIBZ Polyclonal Antibody

Package Size	100 μ g (400 μ L / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20 $^{\circ}$ C Once thawed, store at 4 $^{\circ}$ C. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of mouse CIBZ by peptide affinity chromatography.
Working dilution	For Western blotting ; 3 μ g/mL



Western blotting

Sample : ① Myc-CIBZ protein
② mouse brain lysate
③ mouse kidney lysate

Preparation of antibodies and instruction :

Dr. Kawaichi M. at Division of Gene Function in Animals, Graduate School of Biological Sciences, Nara Institute of Science and Technology



Code No.KAL-KR103

Anti Mouse CIBZ Polyclonal Antibody

【Reference】

1. Sasai N. et al. : Genes Cells. 2005 Sep;10(9):871-85.
2. Kiefer H. et al. : Mol Cell Biol. 2005 Mar;25(5):1713-29
3. Melnick A. et al. : Mol Cell Biol. 2002 Mar;22(6):1804-18
4. Deltour S. et al. : Mol Cell Biol. 2002 Jul;22(13):4890-901

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>


Phone : +81-3-5632-9617

e-mail : export@cosmobio.co.jp

FAX : +81-3-5632-9618

抗 マウス CIBZ ポリクローナル抗体

容量	100 μ g (400 μ L/vial)
形状	ウサギポリクローナル抗体 0.25mg/mL、凍結品
バッファー	PBS [2%ブロッカー(安定化蛋白)、0.1%proclin 含有]
保管方法	−20℃以下 抗体を低濃度にて冷蔵保管されますと、失活する恐れがあります。 融解後は4℃で保存し、お早めにご使用下さい。 また凍結融解を繰り返すことは避けて下さい。
製造方法	マウス CIBZ の部分ペプチドを免疫して得られたウサギの抗血清より、ペプチドアフィニティーカラムにて精製。
使用濃度	ウェスタンブロッティング: 3 μ g/mL



ウェスタンブロッティング

Sample : ① Myc-CIBZ タンパク
 ② マウス脳 lysate
 ③ マウス腎臓 lysate

提供 : 奈良先端科学技術大学院大学
 バイオサイエンス研究科
 動物遺伝子機能学講座
 川市 正史 先生
 松田 永照 先生

抗 マウス CIBZ ポリクローナル抗体

【参考文献】

1. Sasai N. et al. : Genes Cells. 2005 Sep;10(9):871-85.
2. Kiefer H. et al. : Mol Cell Biol. 2005 Mar;25(5):1713-29
3. Melnick A. et al. : Mol Cell Biol. 2002 Mar;22(6):1804-18
4. Deltour S. et al. : Mol Cell Biol. 2002 Jul;22(13):4890-901



人と科学のステキな未来へ

コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
URL : <http://www.cosmobio.co.jp/>

● 営業部（お問い合わせ）

TEL : (03) 5632-9610 FAX : (03) 5632-9619

TEL : (03) 5632-9620