




Code No.KAL-KR111

For research use only

Anti human Tensin2 Polyclonal Antibody

Package Size	50µg (200µL / vial)
Format	Rabbit polyclonal antibody 0.25mg/mL
Buffer	PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage	Store below -20°C Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.
Purification method	This antibody was purified from rabbit serum immunized with synthesized partial peptide of human Tensin2 by peptide affinity chromatography.
Working dilution	For Immunohistochemistry ; 5µg/mL For Western blotting ; 2µg



Immunohistochemistry

Sample : Mouse Kidney

Preparation of antibodies and instruction :

Dr.Agui T. and Sasaki N.

The Graduate School of Veterinary Medicine, Hokkaido University Department of Disease Control

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

Phone : +81-3-5632-9617

e-mail : export@cosmobio.co.jp

FAX : +81-3-5632-9618



Code No.KAL-KR111

Anti human Tensin2 Polyclonal Antibody

【Reference】

1. Lo SH. : Int J Biochem Cell Biol. 2004 Jan;36(1):31-4.
2. Lo SH. et al. : J Cell Biol. 1997 Mar 24;136(6):1349-61.
3. Calderwood DA et al. : Proc Natl Acad Sci U S A. 2003 Mar 4;100(5):2272-7.
4. Chen H. : Proc Natl Acad Sci U S A. 2002 Jan 22;99(2):733-8.
5. Lo SH. : Cancer Res. 2002 Aug 1;62(15):4217-21.

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>


Phone : +81-3-5632-9617

e-mail : export@cosmobio.co.jp

FAX : +81-3-5632-9618

抗 ヒト Tensin2 ポリクローナル抗体

容量	50 μ g (200 μ L/vial)
形状	ウサギポリクローナル抗体 0.25mg/mL、凍結品
バッファー	PBS [2%ブロッカーエース(安定化蛋白)、0.1%proclin 含有]
保管方法	−20℃以下 抗体を低濃度にて冷蔵保管されますと、失活する恐れがあります。 融解後は4℃で保存し、お早めにご使用下さい。 また凍結融解を繰り返すことは避けて下さい。
製造方法	ヒト Tensin2 の部分ペプチドを免疫して得られたウサギの抗血清より、ペプチドアフィニティーカラムにて精製。
使用濃度	免疫組織染色: 5 μg/mL ウェスタンブロッティング: 2 μg/mL



免疫組織染色

Sample : マウス(C57BL6/J) 腎臓

提供 : 北海道大学大学院 獣医学研究科
 安居院 高志 先生
 佐々木 宣哉 先生



人と科学のステキな未来へ

コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
 URL : <http://www.cosmobio.co.jp/>

● 営業部 (お問い合わせ)

TEL : (03) 5632-9610 FAX : (03) 5632-9619

TEL : (03) 5632-9620

抗 ヒト Tensin2 ポリクローナル抗体

【参考文献】

1. Lo SH. : Int J Biochem Cell Biol. 2004 Jan;36(1):31-4.
2. Lo SH. et al. : J Cell Biol. 1997 Mar 24;136(6):1349-61.
3. Calderwood DA et al. : Proc Natl Acad Sci U S A. 2003 Mar 4;100(5):2272-7.
4. Chen H. : Proc Natl Acad Sci U S A. 2002 Jan 22;99(2):733-8.
5. Lo SH. : Cancer Res. 2002 Aug 1;62(15):4217-21.



人と科学のステキな未来へ

コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
URL : <http://www.cosmobio.co.jp/>

● 営業部（お問い合わせ）

TEL : (03) 5632-9610 FAX : (03) 5632-9619

TEL : (03) 5632-9620