

MOUSE ANTI INSLIN-LIKE GROWTH FACTOR-1 [IGF-1] MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-021
Clone	27E10
Concentration	1 mg/ml
Volume	100 UL
Subclass	lgG1k
Immunogen	Recombinant Insulin-like growth factor -1
Specificity	Human Insulin-like growth factor -1
	No cross reactivity with Human Insulin-like growth factor -2
Applications	Sandwich ELISA (This antibody can be used as a capture antibody
	in a human Insulin-like growth factor -1 ELISA in combination with
	biotinlated, monoclonal antibody 16G4.)
Format	Purified immunoglobulin from BALB/C Mouse Ascites
	Anmmonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl,
	0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months.
	Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic



抗ヒト IGF-1 モノクローナル抗体

- 品番 KR-021
- クローン 27E10
- 型 IgG1, ĸ
- 由 来 BALB/Cマウス腹水
- 精 製 法 硫安塩析, r プロテインA
- 溶 液 PBS (0.09%アジ化ナトリウム含む)
- 数量 100µg
- 保存法 4~8℃ 長期保存の場合は、-70℃以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
- 備考抗ヒト IGF-I 抗体(27E10) はサンドイッチ ELISA 法の固相に適しています。なお、二次抗体としては抗ヒト IGF-I 抗体(16G4) が適しています。ヒト IGF-II との交差性はありません。

