

**MOUSE ANTI MYOGLOBIN  
MONOCLONAL ANTIBODY**

<b>Catalog No.</b>	LNM-KR-027
<b>CLONE:</b>	No.29
<b>QUANTITY:</b>	0.1mg
<b>SPECIFICITY:</b>	Human Myoglobin
<b>IMMUNOGEN:</b>	Human Myoglobin
<b>Type/Subtype:</b>	IgG1 $\kappa$
<b>APPLICATIONS:</b>	Sandwich ELISA (This antibody can be used as a capture antibody in a human Myoglobin ELISA in combination with biotinlated, monoclonal antibody No.44.)
<b>FORMAT:</b>	Purified immunoglobulin from BALB/C Mouse Ascites (Affinity chromatography)
<b>PRESENTATION:</b>	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
<b>STORAGE/HANDLING:</b>	2 to 8°C, Avoid repeated freeze/thaw cycle.

***For research use only; not for use as a diagnostic.***



**COSMO BIO CO., LTD.**  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN  
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com