## MOUSE ANTI MYOGLOBIN MONOCLONAL ANTIBODY

Catalog No. LNM-KR-028

CLONE: No.44

**QUANTITY:** 0.1mg

SPECIFICITY: Human Myoglobin

IMMUNOGEN: Human Myoglobin

Type/Subtype:  $IgG1 \kappa$ 

**APPLICATIONS:** Sandwich ELISA (This antibody can be used as a detection

antibody in a human Myoglobin ELISA combination with

the monoclonal capture antibody No.29.)

FORMAT: Purified immunoglobulin from BALB/C Mouse Ascites

(Affinity chromatography)

**PRESENTATION:** Liquid in 15mM Phosphate buffer containing 0.15M NaCl,

0.09% sodium azide, pH 7.2

**STORAGE/HANDLING:** 2 to 8°C, Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic.

