



MOUSE ANTI PEPSINOGEN 2 MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-032
Clone	2D5
Concentration	1 mg/ml
Volume	100 UL
Subclass	IgG1k
Immunogen	Human Pepsinogen 2
Specificity	Human Pepsinogen 2
Applications	Sandwich ELISA (This antibody can be used as a detection antibody in a human Pepsinogen II ELISA combination with the monoclonal capture antibody 10E11.)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Anmmonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic



抗ヒトペプシノゲンⅡモノクローナル抗体

品 番	KR-032
クローン	2D5
型	IgG1, κ
由 来	BALB/Cマウス腹水
精製法	硫酸塩析, r-プロテインA
溶 液	PBS (0.09%アジ化ナトリウム含む)
数 量	100 μ g
保 存 法	4~8℃ 長期保存の場合は、-70℃以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒトペプシノゲンⅡ抗体(2D5)はサンドイッチELISA法の二次抗体に適しています。なお、固相としては抗ヒトペプシノゲンⅡ抗体(10E11)が適しています。



コスモ・バイオ株式会社