

抗ヒトプレアルブミンモノクローナル抗体

品 番	KR-064
クローン	No. 21
型	IgG1, κ
由 来	BALB/Cマウス腹水
精製法	硫酸塩析, r-プロテインA
溶 液	PBS (0.09%アジ化ナトリウム含む)
数 量	100 μ g
保 存 法	4~8°C 凍結融解は出来るだけ避けて下さい。
備 考	抗ヒトプレアルブミン抗体 (No.21) はサンドイッチ ELISA 法に適しています。



人と科学のステキな未来へ

コスモ・バイオ株式会社

— 商品に関するお問い合わせ —

TEL: 03-5632-9610 FAX: 03-5632-9619

本社所在地 〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

**MOUSE ANTI prealbumin
MONOCLONAL ANTIBODY**

Catalog No. LNM-KR-064

CLONE: No.21

QUANTITY: 0.1mg

SPECIFICITY: Human prealbumin

IMMUNOGEN: Human prealbumin

Type/Subtype: IgG1, κ

APPLICATIONS: Sandwich ELISA(This one clone is useful for both labeled antibody and capture antibody.)

FORMAT: Purified immunoglobulin from BALB/C Mouse Ascites
(Affinity chromatography)

PRESENTATION: Liquid in 15mM Phosphate buffer containing 0.15M NaCl,
0.09% sodium azide, pH 7.2

STORAGE/HANDLING: 2 to 8°C, Avoid repeated freeze/thaw cycle.

For research use only; not for use as a diagnostic.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com