



## Anti (chMAP) Myelin protein zero (P0)

<b>Form</b>	Whole serum
<b>Host</b>	Rabbit
<b>Antigen</b>	Multiple Antigen Peptide*
<b>Application notes</b>	ELISA/ IHC/ WB <b>Recommended use and dilutions</b> According to the ELISA assay, results are positive for dilutions up to 200,000 fold against linear synthetic peptides For Western Blotting, dilute > 1 :1,000 For Immunofluorescence, dilute> 1 :500 (1/200 – 1/2000) Optimal diluting should be determined by the user.
<b>Cross Reactivity</b>	HU/ MS/ RAT/ CHK
<b>Volume</b>	100 ul
<b>Storage</b>	Store below -40°C (Store below -80°C for long term storage) After thawing, Avoid freeze-thraw.
<b>Stability</b>	For 3 years at -80C Preservatives such as NaN <sub>3</sub> are not included

\*The cross reactivity is different depending on the species. Please confirm the dilution ratio.

