



COSMO BIO Co., LTD.

Inspiration for Life Science

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. LSL-LB-5199

Anti Mite Extract

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
	This product does not contain preservatives such as NaN3.
Volume	100µl
Concentration	
Specificity	Mite Crude Extract
Antigen Isotype	Mite Crude Extract (#LG-5119)

Application notes	ELISA, WB, IHC, IF Recommended use Positive to 1/25,000 Mite Extract-Df (#LG-5339) and to 1/50,000 Mite Extract-Dp (#LG-5449) and to 1/3,000 Mite Extract-Tp (#LG-5559) . Immunohistochemistry on Frozen or Paraffin section Recommended dilutions Immunofluorescence; More than 1/50 (recommended: 1/50~1/500) Western Blotting; More than 1/500 Optimal dilutions/concentrations should be determined by the end user. Staining Pattern
Cross reactivity	Mite Cross reaction will differ between species.
Storage	Store below -20°C (below -70°C for prolonged storage). After thawing, store in small aliquots in sealable vials and store below -70°C. To prevent degradation from repeated thawing, store the antiserum between 0 to 4°C after second thawing. Stable for three years at -70°C
References	

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Anti Mite Extract.

I. 内容

lot No. 774100

免疫動物	ウサギ
性状・包装サイズ	全血清・100 μ l
抗原由来	Mite Crude Extract (LG-5119)
特徴	Mite Extract-Df (LG-5339) に対して25,000倍希釈、 Dp (LG-5449) に対して50,000倍希釈、Tp (LG-5559) に対して3,000倍希釈まで陽性。

II. 保存上の注意

-20℃以下（長期間の場合は-70℃以下）で凍結して下さい。
 解凍後は密栓のできる小型容器に研究の規模に応じて少量ずつ分注し、
 -70℃以下で保存して下さい。
 凍結融解の繰り返しによる力価の低下を避ける為、再解凍後の抗血清
 は0~4℃に保ち操作・保存して下さい。

III. 安定性

-70℃で3年間安定。
 但し、NaN₃等の防腐剤は入っていません。

製造元

総発売元

 株式会社 エル・エス・エル

 コスモ・バイオ株式会社