

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. LSL-LB-5509

Anti Rhodopsin

Product type Primary antibodies

HostRabbitSourceSerumFormLiquid

This product does not contain preservatives such as NaN3.

Volume 100µl

Concentration

Specificity Rhodopsin **Antigen** Octopus Retina

Isotype

Application notes ELISA, WB, IHC

Recommended use

According to the ELISA assay, results are positive for dilutions up to 500,000 fold.

Immunohistochemistry on Frozen or Paraffin section

Recommended dilutions

Immunofluorescence; More than 1/500 (recommended: 1/200~1/10,000)

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Appears as band corresponding Rhodopsin by Western blotting.

Cross reactivity Cross react with Octopus, Calamary and frog.

Low cross react with chicken and mammalian.

Cross reaction will differ between species.

Storage Store below -20°C (below -70°C for prolonged storage). After thawing, store in small

aliquots in sealable vials and store below -70°C. To prevent degradation from repeated

thawing, store the antiserum between 0 to 4°C after second thawing.

Stable for three years at -70°C

References

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp e-mail: export@cosmobio.co.jp

[Dutside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: +81-3-5632-9610

FAX: +81-3-5632-9618 FAX: +81-3-5632-9619

Anti (octopus) Rhodopsin

I. 内容

lot No. 842041

免 疫	動	物	ウサギ
性状・包装サイス			全血清・100μℓ
カ	\$ 100 K T	価	ELISA で500,000 倍希釈まで陽性
抗原	由	来	ミズダコ・網膜
種間	交	差	<u>タコ</u> ・イカ・カエルと交差*
特		徴	網膜抽出物に対してWestern blottingでRhodopsinに対応するバンドと反応。組織染色で網膜の桿細胞(ロドクシンの存在部位)が染色される。
標準希釈率			蛍光抗体法で1:500以上 (1/200 ~1/10,000)

- *種によって交差性が異なります、希釈率等充分検討下さい。ニワトリ及びホ乳動物との交差性は低いと思われる。
- ()内は推薦希釈倍率

Ⅱ. 保存上の注意

-20℃以下(長期間の場合は-70℃以下)で凍結して下さい。 解凍後は密栓のできる小型容器に研究の規模に応じて少量づつ分注し、-70℃以下で保存して下さい。 凍結融解の繰り返しによる力価の低下を避ける為、再解凍後の抗血清 は0~4℃に保ち操作・保存して下さい。

Ⅲ. 安定性

-70℃で3年間安定。 但し、NaN。等の防腐剤は入っていません。

製造元

総発売元

SI 株式会社 エル・エス・エル

コスモ・バイオ株式会社