



人と科学のステキな未来へ

コスモ・バイオ株式会社

MAB Institute, Inc.
<http://www.monoclo.com>

Anti-RNA Polymerase 2 CTD, mouse monoclonal antibody

Catalog No. MABI0601-100

Lot No. 18028

Product Description

Clone No.	MABI0601
Host	Mouse
Source	Culture supernatant of serum free medium
Isotype	IgG1
Size	100 μ l
Concentration	1 mg/ml
Antigen	18 amino acid residues of human RNA polymerase CTD,
Buffer	PBS+0.05% sodium azide
Storage	Store below -20°C
Expiration Date	Oct 2021

Application notes

Immunoprecipitation	1-5 μ g/5 μ l Sepharose
Immunostaining	0.5-1 μ g/ml
Immunoblot	0.5-1 μ g/ml

Attention

Mouse IgG has weak combination with protein A and protein G. Please use the beads which connect anti-mouse IgG antibodies such as anti-mouse IgG sepharose for immunoprecipitation.



人と科学のステキな未来へ

コスモ・バイオ株式会社

MABI 0601

