



人と科学のステキな未来へ

コスモ・バイオ株式会社

MAB Institute, Inc.

<http://www.monoclo.com>

## Anti-phospho RNA polymerase2 CTD (Ser5), mouse monoclonal antibody

Catalog No. MABI0603-100

Lot No. 18028

### Product Description

Clone No. MABI0603

Host Mouse

Source Culture supernatant of serum free medium

Isotype IgG1

Size 20 $\mu$ l, 100  $\mu$ l

Concentration 1 mg/ml

Antigen 18 amino acid residues of human RNA polymerase CTD, phosphorylated at Ser5

Buffer PBS + 0.05% sodium azide

Storage Store below -20°C

Expiration Date Oct 2021

### Application notes

Immunoprecipitation 1-5  $\mu$ g/5  $\mu$ l Sepharose

Immunostaining 0.5-1  $\mu$ g/ml

Immunoblot 0.5-1  $\mu$ g/ml

### Attention

Mouse IgG has weak combination with protein A and protein G. Please use the beads which connect anti-mouse IgG antibodies such as anti-mouse IgG sepharose for immunoprecipitation.



人と科学のステキな未来へ

# コスモ・バイオ株式会社

MABI 0603

