



## Monocyte Wash Medium

### Product Description

Monocyte Wash Medium (PMC-BMMW-COS) is formulated to allow for high viability and optimal growth of monocytes after thawing and upon plating. It is recommended for use with Monocyte Culture Medium (Code No. PMC-BMMG-COS). This product is an optional part of Monocyte Culture Kit (Rat, PMC-BMM01-COS).

### Components / Storage

Product name / Code No.	Size	Quantity	Storage Conditions	Shelf Life
Monocyte Wash Medium (Code No. PMC-BMMW-COS)	50 mL	1	-20°C Freezer	6 months
			4°C	3 months

# Shipping: dry ice

### Precautions

This product is for research use only, not for human or diagnostic use

### References

- (1) Sunao Takesita, Keisuke Kaji, Akira Kudo. Identification and Characterization of the New Osteoclast Progenitor with Macrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. JOURNAL OF BONE AND MINERAL RESEARCH Volume 15, Number 8(2000) .1477-1488.

Manufactured by **Primary Cell Co., Ltd.**



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619