

# Differentiation medium for Chondrocyte Culture kit (Rabbit)

# (500mL, Code No. PMC-CHCM-COS)

#### For research use only

### **Product Description**

Differentiation medium for Chondrocyte Culture kit (Rabbit) (Code No. PMC- CHCM04-COS) is complete medium optimized for *in vitro* differentiation of cells isolated from rabbit articular cartilage into mature chondrocyte. This product is available for optional Chondrocyte Culture Kit (Rabbit) (PMC-CHCM-COS).

### **Components/ Storage**

Differentiation Medium for Chondrocyte Culture kit (Rabbit) (Code No. PMC- CHC04-COS) is sterile, liquid basal medium (RPMI1640) which contain essential and non-essential amino acids, vitamins, other organic compounds, trace minerals, inorganic salts, FBS, and antibiotics.

Product Name/Code No.	Size	Quantity	Storage Conditions	Stability
Differentiaiton Medium for Chondrocyte	500 mL	1	-20°C Freezer	6 months
Culture kit (Rabbit)			4°C	3 months
(Code No. PMC- CHCM-COS)				

\* Shipping: dry ice

## Precautions

This product is for research use only, not for human or diagnostic use.

