



COSMO BIO CO., LTD.
Inspiration for Life Science

[Cat. No. MPCDNA-DT-COS]

Mouse Pancreas First Strand cDNA (dT)

For research use only

Product information

Lot Number: ***

Source: ICR Mouse Pancreas

Synthesis format: Mouse Pancreas cDNA is the first strand cDNA synthesized from total RNA 0.5µg using oligo (dT)₂₀ primer.

Volume: 50µL

Reagent: 2mM Tris-HCl (pH 8.4), 5mM KCl, 0.5mM MgCl₂, 1mM DTT

Application: This product can be used for RT-PCR. 0.125µL cDNA is available for detecting high expression gene (Insulin, ribosomal RNA). 2.5µL cDNA is available for detecting normal expression gene (β-actin, Pdx-1).

Storage: ≤-20°C

shipment: dry ice

Precautions

- Read the instructions carefully before using.
- This product is for research use only, not for human or diagnostic use.
- Always wear gloves and lab coat when handling.

Quality Control

RNA Integrity Number (RIN): *** (detected by 2100 Bioanalyzer, Agilent Technologies, Inc.)

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN
TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618 / e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com