

Anti human PXR common mouse monoclonal antibody

PXR: Pregnane X receptor, SXR

Code No PP-H4417-00

Clone No. H4417

Lot. ***

Concentration 1 mg/mL

Volume 100 uL

Ig Class G2a

Description Pregnane-activated receptor (PXR,SXR, PAR, PAR1, PAR2, NR112) is a member of nuclear receptor subfamily , which each orthologues historically given different names as pregnane-activated receptor in mice and steroid- and xenobiotic-sensing nuclear receptor in human, PXR binds to rifampicin (an antibiotics) is the most efficient activator in human. Various studies revealed PXR regulates CYP3A gene expression as well as other xenobiotic metabolisms, such as oxidation, conjugation and transport. Many chemicals are known to bind for PXR as activators, eg. the HIV protease inhibitor ritonavir, the anticancer drug paclitaxel, the endocrine disruptor bisphenol A. Expression of PXR founds in the liver, small intestine and colon in the human, rabbit and mouse where CYP3A genes are expressed or induced.

Nomenclature NR11 2

Genbank AF084645

Origin Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human PXR-1(1-40 aa) .

Specificity This antibody specifically recognizes human PXR-1 and PXR-2. Not yet tested in other species.

Purification Ammonium sulfate fractionation

Formulation Physiological saline with 0.1% NaN₃ as a preservative.

Application / Recommended Concentration

In order to obtain the best results, optimal working dilutions should be determined by each individual user.

Western Blot 1ug/mL

Non reducing Western Blot 3ug/mL

ELISA 2ug/mL (A450=1.0)

Immunoprecipitation Decide by use

Supershift Assay Not yet tested

Chromatin immunoprecipitation Not yet tested

Immunohistochemistry Not yet tested

Storage Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.

Reference

Notes Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited.

MADE IN JAPAN

Aug 29, 2006