

Anti human PPAR delta mouse monoclonal antibody

PPAR delta: Peroxisome Proliferator-Activated Receptor delta

Code No	PP-K9436-00	Application / Recommended Concentration In order to obtain the best results, optimal working dilutions should be determined by each individual user.	
Clone No.	K9436	Western Blot	1 ug/mL
Lot.	***	Non reducing Western Blot	Not yet tested
Concentration	1 mg/mL		
Volume	100 uL	ELISA	0.5 ug/mL
lg Class	G2a	Immunoprecipitation	Not yet tested
Description	Peroxisome proliferator-activated receptor beta (PPARb, NUC1, PPARd; NR1C2) is a member of orphan nuclear receptor. PPARb is expressed in uterus and placenta. It plays important roles in lipid and glucose metabolism, and have been implicated in obesity-related metabolic diseases such as hyperlipidemia, insulin resistance, and coronary artery disease. Three members were called PPARa, b, g. RXR is an obligate partner for PPAR.	Supershift Assay	100 ug/mL
		Chromatin immunoprecipitation	Not yet tested
		Immunohistochemistry	Not yet tested

Nomenclature	NR1C2			
Genbank	L07592			
Origin	Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a BALB/c mouse immunized with Baculovirus-expressed recombinant human PPAR delta (1-60 aa).	Storage	Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.	
Specificity	This antibody specifically recognizes human PPAR delta but does not recognize human PPAR alpha and gamma. Not yet tested in other species.	Reference		
Purification	Ammonium sulfate fractionation			
Formulation	Physiological saline with 0.1% NaN3 as a preservative.	Notes	Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.	

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited. **MADE IN JAPAN**

Aug 23, 2006



Distributor:

2-20,Toyo, 2Chome,Tokyo 135-0016,Japan TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 http://www.cosmobio.co.jp



Manufactured by Perseus Proteomics Inc. 4-7-6, Komaba, Meguro-ku, Tokyo 153-0041, Japan http://www.ppmx.com