

Anti human PTX3 mouse monoclonal antibody

Pentraxin 3

Code No	PP-PPZ1272-00
Clone No.	PPZ1272
Lot.	A-2
Concentration	1 mg/mL
Volume	100 uL
Ig Class	G1
Description	Pentraxins are a superfamily of conserved proteins characterized by the pentraxin domain. CRP and SAP are recognized as classical short pentraxins, whereas Pentraxin3(PTX3) belongs to the long pentraxins. CRP and SAP are induced in the liver in response to IL-6. In contrast, PTX3 is produced by a variety of tissues and cells, such as vascular endothelial cells, macrophages, and neutrophils, predominantly in response to proinflammatory signals (bacterial products, IL-1, and TNF).
Genbank	BC039733
Origin	Produced in BALB/c mouse ascites after inoculation with hybridoma of mouse myeloma cells (NS-1) and spleen cells derived from a mouse immunized with recombinant human PTX3 (18-381 aa).
Epitope	80-108 a.a.
Specificity	This antibody specifically recognizes human PTX3. Not yet tested in other species.
Purification	Ammonium sulfate fractionation.
Formulation	Physiological saline with 0.1% NaN ₃ as a preservative.

Application / Recommended Concentration

In order to obtain the best results, optimal working dilutions should be determined by each individual user.

Western Blot 1 ug/mL

Non reducing Western Blot 1 ug/mL

ELISA Decide by use

Immunoprecipitation Not yet tested

Supershift Assay Not yet tested

Chromatin immunoprecipitation Not yet tested

Immunohistochemistry Not yet tested

Storage

Store at 2 - 8 °C up to one month. For long-term storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is not recommended. Storage in a frost-free freezer is not recommended.

Reference

Notes

Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

Not for Diagnostic or Therapeutic use. Purchase of this product does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written consent of Perseus Proteomics Inc. is prohibited.

MADE IN JAPAN

Nov 18, 2008



COSMO BIO CO., LTD.
Inspiration for Life Science

Distributor:
Cosmo Bio Co.,Ltd.
2-20,Toyo, 2Chome,Tokyo 135-0016,Japan
TEL:+81-3-5632-9617 FAX:+81-3-5632-9618
<http://www.cosmobio.co.jp>



Manufactured by
Perseus Proteomics Inc.
4-7-6, Komaba, Meguro-ku, Tokyo 153-0041
Japan
<http://www.ppmx.com>

Anti human PTX3 mouse monoclonal antibody

Pentraxin 3

製品コード PP-PPZ1272-00

Clone No. PPZ1272

Lot. A-2

濃度 1mg/mL

容量 100 μ L

Ig class G1

Genbank BC039733

由来 リコンビナントヒトPTX3(18-381 a.a.)を免疫したマウスの脾臓細胞と、マウスミエローマ細胞(NS-1)を融合して得たハイブリドーマをBALB/c マウスに接種して得られた腹水

Epitope 80-108 a.a.

特異性 ヒトPTX3に特異的に反応する。その他の動物種との交差反応は未検討。

精製法 硫酸塩析法

溶媒 生理的食塩水(防腐剤として0.1% NaN3添加)

Application 使用濃度は実験にあわせて至適化が必要です。

Western Blot 可
参考使用濃度 1 μ g/mL

非還元 Western Blot 可
参考使用濃度 1 μ g/mL

ELISA 可
参考使用濃度 適宜調整してください

免疫沈降 未検討
参考使用濃度 -

Supershift Assay 未検討
参考使用濃度 -

クロマチン免疫沈降 未検討
参考使用濃度 -

免疫染色 未検討
参考使用濃度 -

保存方法 1ヶ月程度の保存の場合は、2~8 $^{\circ}$ Cで保存可能です。長期保存の場合は、抗体を小分けした上で、-20 $^{\circ}$ C以下での保存をお勧めします。また、凍結融解を繰り返すと、抗体が劣化し、本来の性能が得られない場合があるため、お避けください。

参考文献

備考 溶媒に含まれるNaN₃は、鉛や銅と反応し爆発性化合物を形成する恐れがあります。廃棄の際には大量の水と一緒に希釈して廃棄してください。

FOR RESEARCH ONLY. NOT FOR USE IN HUMANS.

本製品は研究目的のみで使用することができます。人やそのほか動物の疾病診断、治療・予防に使用することはできません。研究目的以外で使用する場合は、あらかじめ弊社にご相談ください。

MADE IN JAPAN

Nov 26, 2008



COSMO BIO CO., LTD.
Inspiration for Life Science

【販売元】
コスモ・バイオ株式会社
〒135-0016 東京都江東区東陽2-2-20
TEL : 03 (5632) -9610 FAX : 03 (5632) -9619
<http://www.cosmobio.co.jp>



【製造】
株式会社 ペルセウスプロテオミクス
〒153-0041 東京都目黒区駒場4-7-6
<http://www.ppmx.com>