



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Human ZNF589 Polyclonal Antibody, Rabbit

Code # : KB4343GNP [FSR02132], (Lot# XXX)

Gene : ZNF589 (Zinc finger protein 589, Stem cell zinc finger protein 1)

Immunogen : GX5380 (GST-fusion protein, 147 amino acids)
GSQQFLSQDELHNHPIPGFHAGNQLHPGNPCPEDQPQSQHPSDKNHRGAEAEQDRVEGGVR
PLFWSTNERGALVGFSSLFQRPPISWGGNRILEIQLSPAQNASSEEVDRISKRAETPGFGAVT
FGECALAFNQKSNLFRQKAVTA

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : This antibody detects human ZNF589 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>