



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Human TRIM28 Polyclonal Antibody, Rabbit

Code # : KB7016GNP [FSR01573], (Lot# XXX)

Gene : TRIM28 (Transcription intermediary factor 1-beta, TIF1-beta, Tripartite motif-containing protein 28, Nuclear corepressor KAP-1, KRAB-associated protein 1, KAP-1, KRAB-interacting protein 1, KRIP-1, RING finger protein 96)

Immunogen : GX5066 (GST-fusion protein, 166 amino acids)
TKSAEAFGKIVAERPSTNSTGPAPMAPPAPGPLSKQSGSSQPMEVQEGYGFSGDDPYSS
AEPHVSGVKRSRSGEGEVSGLMRKVPRVSLERLDLDTADSQPPVFKVFPGSTTEDYNLIVIER
GAAAAATGQPGTAPAGTPGAPPLAGMAIVKEEETEAIGA

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : This antibody detects human TRIM28 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>