



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Human ANP32A Polyclonal Antibody, Rabbit

**Code # :** KB7264GNP [FSR02366], (Lot# XXX)

**Gene :** ANP32A (Acidic leucine-rich nuclear phosphoprotein 32 family member A, Potent heat-stable protein phosphatase 2A inhibitor I1PP2A, Acidic nuclear phosphoprotein pp32, Leucine-rich acidic nuclear protein, Lanp, Putative HLA-DR-associated protein I, PHAPI, Mapmodulin)

**Immunogen :** GX5388 (GST-fusion protein, 153 amino acids)  
LELRNRTPSDVKELVLDNSRSNEGKLEGLTDEFEELEFLSTINVGLTSIANLPKLNKLNKLELSDN  
RVSGGLEVLAEKCPNLTHLNLSGNKIKDLSTIEPLKKLENLKSLLDLFNCEVTNLNDYRENVFKLLP  
QLTYLDGYDRDDKEAPDSDAE

**Format :** Rabbit Whole Serum

**Preservative :** 0.05 % NaN<sub>3</sub>

**Quantity :** 100 µL

**Specificity :** This antibody detects human ANP32A protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at -20°C. Avoid freeze-thaw cycles.

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**  
2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>