



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Human ZFAND5 Polyclonal Antibody, Rabbit

Code # : KB7270GNP [FSR02370], (Lot# XXX)

Gene : ZFAND5 (Acidic leucine-rich nuclear phosphoprotein 32 family member A, Potent heat-stable protein phosphatase 2A inhibitor I1PP2A, Acidic nuclear phosphoprotein pp32, Leucine-rich acidic nuclear protein, Lanp, Putative HLA-DR-associated protein I, PHAPI, Mapmodulin)

Immunogen : GX5436 (GST-fusion protein, 173 amino acids)
EHLQRQQNSGRMSPMGTASGSNSPTSDSASVQRADTSLNCEGAAGSTSEKSRNVPVAALP
VTQQMTEMSISREDKITTPKTEVSEPVVTQPSPSVSQPSTSQSEEKAPELPKPKKNRCFMCRK
KVGLTGFDCRCGNLFCGLHRYSDKHNCOPYDYKAEAAAKIRKENPVVVAE

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : This antibody detects human ZFAND5 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>