



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Human NCOR2 Polyclonal Antibody, Rabbit

Code # : KD0137GNPAF [FAF07721], (Lot# XXX)

Gene : NCOR2 (Nuclear receptor corepressor 2, N-CoR2, Silencing mediator of retinoic acid and thyroid hormone receptor, SMRTe, SMRT, Thyroid-, retinoic-acid-receptor-associated corepressor, T3 receptor-associating factor, TRAC, CTG repeat protein 26, SMAP270)

Immunogen : GX5151 (GST-fusion protein, 169 amino acids)
GSTQPVAQTWRAATEPRYPHSLSYVQIARTHTDVGLEYYQHHSRDYASHLSPGSIQPPQRRR
PSLLSEFQPGNERSQELHLRPESHSYLPELGKSEMEFIESKRPRLELLPDPLLRPSPLLATGQPA
GSEDLTKDRSLTGKLEPVSPSPPHTDPELELVPPRLSKEE

Format : Affinity Purified Rabbit IgG

Constitution : PBS containing with 40% glycerol and 0.02% of NaN₃

Quantity : 50 µg (250 µL)

Specificity : This antibody detects human NCOR2 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>