

Anti-Mouse Rasal3 (FLJ00087) Polyclonal Antibody, Rabbit

Code # : MFL0087 [mFLJ00087], (Lot# 001)

Gene : mouse Rasal3 (RAS protein activator like 3, mFLJ00087)

Immunogen : GX1975 (156 amino acids)
LPRKPSVPWQRQMDQPGDRYQTTGTHRPVVGKLAIEIQCEVAIFREAQKALSLLVESLSTQV
QALKEQQEHFRCQLQDLYSRLGAGISKLDSKGGLPSNGSHRLKSLEQRLTEMECSQDQLR
DSLQSLQLLSKTPGSRSQPLPLKAPCVNGADLSMGT

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX1975. This antibody detects mouse Rasal3 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

