



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse PNPLA7 (FLJ00415), Rabbit Polyclonal

Code # : MFL0415 (Lot# xxx)

Gene : patatin-like phospholipase domain containing 7 (PNPLA7), transcript variant 1, mRNA

Immunogen : GX2403 (238 amino acids)
ARSMGAKVVIAIDVGSRDETLTNYGDALSGWLLWKRWNPLATKVKVLNMAEIQTRLAY
VCCVRQLEMVKNSDYCEYL RPPIDSYRTLDFGKDFEICEVGYQHGRVFDIWRVRSGLVLEK
MLQDQQGTSKRKDCGVFTCPNSSFTDLAEIVSRIEPAKVAAVDDSDYQTEYEEELPAIP
KETYADFQSTGIELSDSEYEPSMLQGPPSLTSPEQSQDSFPWLPNQDDQGPRLEHPS

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX2403. This antibody detects mouse PNPLA7 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute

2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>