



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Mouse ZSCAN12 (KIAA0426) Polyclonal Antibody, Rabbit

- Code # :** MK04260910, (Lot# XXX)
- Gene :** mouse zinc finger and SCAN domain containing 12 (ZSCAN12)  
(mZSCAN12, mKIAA0426)
- Immunogen :** GX0926 (GST-fusion protein, 238 amino acids)  
PGRKVHGCDECGKSFTQHSRLIEHKRVHTGDRPYKCEVCGKTFWRVRTVLIRHKVVHTGEKPYK  
CNECGRAFGQWSALNQHQRLHSGEKHYHCNECGKAFKAGLFFHHLKSHRRNRPYQCLQC  
NKSFNRRSTLSQHQGVHTGAKPYECNDCGKAFVYNSSLATHQETHHKEKPFQSGPIQQQRN  
HTKEKPYKCSVCGKAFIQKISLIEHEQIHTGERPYKCAEGGKAFIQMSELTEH
- Format :** Affinity Purified Rabbit IgG
- Constitution :** PBS containing with 50% glycerol and 0.02% of NaN<sub>3</sub>
- Quantity :** 50 µg (250 µL)
- Specificity :** Specific to recombinant protein GX0926. This antibody detects mZSCAN12 protein. Details are shown in the InGaP database.  
(<http://webcreate.kazusa.or.jp/create/index.jsp>)  
Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Other applications have not been tested.
- Storage :** Store at -20°C. Avoid freeze-thaw cycles.
- References :** 1) Okazaki, N. *et al.*: *DNA Res.*, 10(1), 35 (2003).

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**  
2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>