



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse CNKSR2 (KIAA0902) Polyclonal Antibody, Rabbit

- Code # :** MK09020910, (Lot# XXX)
- Gene :** mouse Connector enhancer of kinase suppressor of Ras 2 (CNKSR2)
(mCNKSR2, mKIAA0902)
- Immunogen :** GX0244 (GST-fusion protein, 171 amino acids)
VSACDPQDDIQPPEVEEEEEEEEEAAAGENVGEKNENREEKLGDSLQDLYRALEEASLSPLGE
HRISTKMEYKLSFIKRCNDPVMNEKLHRLRLKSTLKAREGEVAIIDKVLDPDLTSKEFQQWKQ
MYLDLFLDICQSTTSNDPLSISSEVDVLTSSLTHTHSYIETHV
- Format :** Affinity Purified Rabbit IgG
- Constitution :** PBS containing with 50% glycerol and 0.02% of NaN₃
- Quantity :** 50 µg (250 µL)
- Specificity :** Specific to recombinant protein GX0244. This antibody detects mCNKSR2 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>)
Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Other applications have not been tested.
- Storage :** Store at -20°C. Avoid freeze-thaw cycles.
- References :** 1) Okazaki, N. *et al.*: *DNA Res.*, 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>