



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse SMARCAD1 (KIAA1122) Polyclonal Antibody, Rabbit

- Code # :** MK11220801 (Lot# XXX)
- Gene :** mouse SWI/SNF-related, matrix-associated actin-dependent regulator of chromatin, subfamily A, containing DEAD/H box 1 (mSMARCAD1, mKIAA1122)
- Immunogen :** GX0323 (GST-fusion protein, 175 amino acids)
DSGKFRALGCILSELKQKGDREVLFQFTMMLDILEVLLKHHQHRYLRLDGKTQISERIHLLIDEFN
TDMDFVFLSTKAGGLGINLTSANVVILHDIDCNPYNDKQAEDRCHRVGQTKVELVIKLISQGTIE
ESMLKINQQKLEQDMTTVDEADEGSMPADIATLLKTSMGL
- Format :** Affinity Purified Rabbit IgG
- Constitution :** PBS containing with 50% glycerol and 0.02% of NaN₃
- Quantity :** 50 µg (250 µL)
- Specificity :** Specific to recombinant protein GX0323. This antibody detects mSMARCAD1 protein. It also recognizes human SMARCAD1 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>)
Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Immunocytochemistry (1 : 1,000), Immunoprecipitation, Chromosomal Immunoprecipitation (ChIP). Other applications have not been tested.
- Storage :** Store at -20°C. Avoid freeze-thaw cycles.
- References :** 1) Hara, Y. et al.: *DNA Res.*, 10, 129 (2003).
2) Koga, H. et al.: *DNA Res.*, 11, 293 (2004).
3) Okazaki, N. et al.: *J. Mol. Biol.*, 382, 257 (2008).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>