



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse E2F3 (KIAA0075) Polyclonal Antibody, Rabbit

- Code # :** MKA0075AF, (Lot# XXX)
- Gene :** mouse Transcription factor E2F3 (mE2F3, mKIAA0075)
- Immunogen :** GX1865 (GST-fusion protein, 165 amino acids)
ENQRLAYVITYQDIRKISGLKDQTVIVVKAPPETRLEVPSIESLQIHLASTQGPIEVYLCPEETETH
RPMKTNNQDHNIPKPTSKDLASNNSGHSDCSVSTANLSPLASPNLLQQTEDQIPSNLEGP
FVNLLPPLLQEDYLLSLGEEEGISDLFDAYDLEKL
- Format :** Affinity Purified Rabbit IgG
- Constitution :** PBS containing with 50% glycerol and 0.02% of NaN₃
- Quantity :** 50 µg (250 µL)
- Specificity :** Specific to recombinant protein GX1865. This antibody detects mE2F3 protein. Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Other applications have not been tested.
- Storage :** Store at -20°C. Avoid freeze-thaw cycles.
- References :** 1) Okazaki, N. *et al.*: *DNA Res.*, 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>