



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Mouse ARHGEF7 (KIAA0142) Polyclonal Antibody, Rabbit

**Code # :** MKA0142PA [FIG00426], (Lot# XXX)

**Gene :** mouse Rho guanine nucleotide exchange factor (GEF) 7  
(mARHGEF7, mKIAA0142)

**Immunogen :** GX0938 (GST-fusion protein, 131 amino acids)  
LSKSPKTMKLLPKRKPERKPSDEEFAVRKSTAALEEDAQILKVIEAYCTSAKTRQTLNSTWQG  
TDLMHNHVLADDDQSSLDLGRSSLSRLEPSDLSSEYDSIWTAAHSYRMGSASRKSCCSYI  
SHQN

**Format :** Rabbit IgG purified with Protein A affinity chromatography.

**Constitution :** PBS containing with 50% glycerol and 0.02% of  $\text{NaN}_3$

**Quantity :** 100  $\mu\text{g}$  (200  $\mu\text{L}$ )

**Specificity :** Specific to recombinant protein GX0938. This antibody detects mARHGEF7 protein. It also recognizes human ARHGEF7 protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at  $-20^\circ\text{C}$ . Avoid freeze-thaw cycles.

**References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**  
2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>