



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Mouse SMG7 (KIAA0250) Polyclonal Antibody, Rabbit

**Code # :** MKA0250 [FSR00477], (Lot# XXX)

**Gene :** mouse Smg-7 homolog, nonsense mediated mRNA decay factor (C. elegans) (mSMG7, mKIAA0250)

**Immunogen :** GX1001 (GST-fusion protein, 178 amino acids)  
TSSKAELNPSVASQETSLYSLFEGTPWSPSLPASSDHSTPASQSPHSSNPSSLPSSPPTHNHN  
SAPFSNFGPIGTPDNRDRPADRWKTDKPMGGFGVDYLSATSSSESSWHQASTPSGTWTG  
HGPSMEDSSAVLMESLKSIWSSSMHPGPSALEQLLMQKQKQQRGQGAMNPPH

**Format :** Rabbit Whole Serum

**Preservative :** 0.05 % NaN<sub>3</sub>

**Quantity :** 100 µL

**Specificity :** Specific to recombinant protein GX1001. This antibody detects mSMG7 protein. It also recognizes human SMG7 protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at -20°C. Avoid freeze-thaw cycles.

**References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**

2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>