



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti- Mouse RALGPS1 (KIAA0351) Polyclonal Antibody, Rabbit

Code # : MKA0351 [FSR00057], (Lot# XXX)

Gene : mouse Ral GEF with PH domain and SH3 binding motif 1
(mRALGPS1, mKIAA0351)

Immunogen : GX0800 (GST-fusion protein, 299 amino acids)
DKTTFEKLDYLSKEDNYKRTRDYIRSLKMVPSIPYLGIIYLLDIYIDSAYPASGSIMENEQRSNQ
MNNILRIIADLQVSCSYDHLTTLPHVQKYLKSVRYIEELQKFVEDDNYKLSLRIEPGSSSPRLVSS
KEDLAGPSAGSSSARFRRRPTCPDTSVAGSLPTPPVPRHRKSHSLGNKGRLYATLGPNWVRVP
VRNSPRTRSCVYSPTSPCTCTVGSSATVPTMEGPLRRKTLLEGRKPALSSWTRYWVVLGAT
LLYYGAKSLRGTDRKHVKSTPGKKVSIVGWMVQLPDDPEHPD

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX0800. This antibody detects mRALGPS1 protein. It also recognizes human RALGPS1 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>