



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse ZBTB43 (KIAA0414) Polyclonal Antibody, Rabbit

Code # : MKA0414, (Lot# XXX)

Gene : mouse zinc finger and BTB domain containing 43 (mZBTB43, mKIAA0414)

Immunogen : GX1879 (GST-fusion protein, 161 amino acids)
EQAWDGMVDVHAAAYDEHQVTESVNTMQTDHSAQPSGAEEEEFQIVEKKVEVEFDEQAEGSSYD
EQVDFYGGSSMEEFSGEKLGGNLIGHKQEAALAAGYSENIEMAMGIKEEASHLGFSATDKLYPC
QCGKSFTHKSQDRHRMSMHLGLRPYGC SVCGKKFKMK

Format : Rabbit Whole Serum.

Constitution : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX1879. This antibody detects mZBTB43 protein. It also recognizes human ZBTB43 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619