



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti- Mouse ARHGEF9 (KIAA0424) Polyclonal Antibody, Rabbit

Code # : MKA0424PA [FIG00073], (Lot# XXX)

Gene : mouse Cdc42 guanine nucleotide exchange factor (GEF) 9
(mARHGEF9, mKIAA0424)

Immunogen : GX0371 (GST-fusion protein, 176 amino acids)
FLFDHQMVLCKKDLIRRDILYYKGRIDMDKYEVIDIEDGRDDDFNVSMKNAFKLHNKETEEVHLF
FAKKLEEKIRWLRAFREERKMVQEDEKIGFEISENQKQAAMTVRKASKQKGVNSARSVPPSY
PPPQDPLNQGQYLVPDGIASQVFEFTEPKRSQSPFWQNFSLTPFKK

Format : Rabbit IgG purified with Protein A affinity chromatography.

Constitution : PBS containing with 50% glycerol and 0.02% of NaN_3

Quantity : 100 μg (200 μL)

Specificity : Specific to recombinant protein GX0371. This antibody detects mARHGEF9 protein. It also recognizes human ARHGEF9 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C . Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>