



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse KIAA0556 Polyclonal Antibody, Rabbit

- Code # :** MKA0556PA [FIG01093], (Lot# XXX)
- Gene :** mouse KIAA0556 (mKIAA0556)
- Immunogen :** GX0407 (GST-fusion protein, 188 amino acids)
LWNYNKSPEDTYRGGVKIAHVSLDGLCVSPAEGFLIRKGPNGCHFDFAEILFGDYLQTRLPPAP
TRRLDAKSLERASMDYEAPLMPCGFIFQFQLLSSWGDPPYYIGLTGLELYDEHGERIPLSQNNIAA
FPDSVNALEGVCGDVRTPDKLIDQVNDTSDGRHMWLAIPILPGLVGTWPWGWTGTCTRWAQ
- Format :** Rabbit IgG purified with Protein A affinity chromatography.
- Constitution :** PBS containing with 50% glycerol and 0.02% of NaN₃
- Quantity :** 100 µg (200 µL)
- Specificity :** Specific to recombinant protein GX0407. This antibody detects mKIAA0556 protein. It also recognizes human KIAA0556 protein. Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Other applications have not been tested.
- Storage :** Store at -20°C. Avoid freeze-thaw cycles.
- References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>