



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse HDAC5 (KIAA0600) Polyclonal Antibody, Rabbit

Code # : MKA0600AF, (Lot# XXX)

Gene : mouse histone deacetylase 5 (HDAC5) (mHDAC5, mKIAA0600)

Immunogen : GX0181 (GST-fusion protein, 143 amino acids)
ARCFGHLTRQLMTLAGGRVVLALLEGGHDLTAICDASEACVSALLSVELQPLDEAVLQQKPSVNA
VATLEKVIIEIQSKHWSCVQRFAAGLGC SLREAQTGEKEEAETVSAMALLSVGAEQAQAVATQE
HSPRPAEEPMEQEPAL

Format : Affinity Purified Rabbit IgG

Constitution : PBS containing with 50% glycerol and 0.02% of NaN_3

Quantity : 50 μg (250 μL)

Specificity : Specific to recombinant protein GX0181. This antibody detects mHDAC5 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C . Avoid freeze-thaw cycles.

References : 1) Okazaki, N. *et al.*: *DNA Res.*, 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>