



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse GIGYF2 (KIAA0642) Polyclonal Antibody, Rabbit

Code # : MKA0642 [FSR01599], (Lot# XXX)

Gene : mouse PERQ amino acid-rich with GYF domain-containing protein 2, GRB10-interacting GYF protein 2, Trinucleotide repeat-containing gene 15 protein) (mGIGYF2, mKIAA0642)

Immunogen : GX0949 (GST-fusion protein, 325 amino acids)
NTATCQSQATLSLAEIQKLEERERQLREEQRRQQRELMKALQQQQQQQQKLSGWSNVN
KPAGTTKSLLEIQQEEARQMKGQQQQQQQQQHQQSNRARNSTHSLNHTSLGNSVWGSIN
TGPSNQWASELVSSIWSNADTKNSNMGFWDDAVKEVGPRNSTNKNKNNASLSKSVGVSNRQ
NKKVEEEEKLLKLFQGVNKAQDGFTQWCEQMLHALNTANNLDVPTFVSFLKEVESPYEVHDYT
RAYLGDTSEAKEFAKQFLERRAKQKVNQQRQQQQQQQDQSVWGMNHSTLHVSFQTNQS
NNQQSNFEAVQSGKKKKKKQ

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX0949. This antibody detects mGIGYF2 protein. It also recognizes human GIGYF2 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>