



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse RAB11FIP3 (KIAA0665) Polyclonal Antibody, Rabbit

Code # : MKA0665 [FSR00598], (Lot# XXX)

Gene : mouse RAB11 family interacting protein 3 (class II) (mRAB11FIP3, mKIAA0665)

Immunogen : GX0044 (GST-fusion protein, 196 amino acids)
CTPCLKANIERLEEEKQKMLDEIEELTQRLSEEQENKRKMGDRLSHERHQFQRDKEATQELIED
LRKQLEHLQLLRLEVEQRRGRSSSLGLQEYNSRARESELEQEVRRRLKQDNRNLKEQNDLNG
QIITLSIQGAKSLFSTSFSESLAAEISSVSRDELMEAIQKQEEINFRLQDYIDRIIVAILETNPSILEVK

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX0044. This antibody detects mRAB11FIP3 protein. It also recognizes human RAB11FIP3 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>