



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Mouse ACIN1 (KIAA0670) Polyclonal Antibody, Rabbit

- Code # :** MKA0670AF, (Lot# XXX)
- Gene :** mouse Apoptotic chromatin condensation inducer in the nucleus (ACIN1) (mACIN1, mKIAA0670)
- Immunogen :** GX0829 (GST-fusion protein, 105 amino acids)  
FKRKISVVSATKGVQAGNSDTEGGQPGRKRRWGASTAATQKKPSITTESLKEAVV  
DLHADDSRISEDETERNGDDGTHDKGLKICRTVTQVVPAEGQENGQRE
- Format :** Affinity Purified Rabbit IgG
- Constitution :** PBS containing with 50% glycerol and 0.02% of NaN<sub>3</sub>
- Quantity :** 50 µg (250 µL)
- Specificity :** Specific to recombinant protein GX0829. This antibody detects mACIN1 protein. Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Other applications have not been tested.
- Storage :** Store at -20°C. Avoid freeze-thaw cycles.
- References :** 1) Okazaki, N. *et al.*: *DNA Res.*, 10(1), 35 (2003).

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**  
2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>