



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse JMJD2A (KIAA0677) Polyclonal Antibody, Rabbit

Code # : MKA0677 [FSR00602], (Lot# XXX)

Gene : mouse jumonji domain containing 2A (mJMJD2A, mKIAA0677)

Immunogen : GX2138 (GST-fusion protein, 298 amino acids)
SANALEEDCCLCSLRGGALQRANDDRWVHVSCAVALIARFVNIAERSPVDVSKIPLPRFKLKC
VFCKRRRKRNAGCCVQCSHGRCPTAFHVSCAQAAGVMMQPDDWPFVVFITCFRHKIPNLERA
KGALLSITAGQKVISKHKNRFRYQCEVVRLTTETFYEVNFDDGFSFDNLYPEDIVSQDCLQLGP
PAEGEVVQVRWTDGQVYGAKFVASHPIQMYQVEFEDGSQVLVVKRDDVYTLDEELPKRVKSRL
SVASDMRFNEIFTEKEVKQEKKRQRVINSRYREDYIEPALYRAIME

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX2138. This antibody detects mJMJD2A protein. It also recognizes human JMJD2A protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute

2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>