



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse TBC1D1 (KIAA1108), Polyclonal Antibody, Rabbit

Code # : MKA1108,(Lot# XXX)

Gene : Tbc 1 protein, Tbc1 domain family, member 1 GTPase activator protein

Immunogen : GX0168 (GST-fusion protein, 174 amino acids)
SEVIFKVALSLLGSHKPLILQHENLETIVDFIKNTLPNLGLVQMEKTISQVFEMDIAKQL
QAYEVEYHVLQEELIESSPLSDNQRMEKLEKTNSSLRKQNLDLLLEQLQVANARIQSLEAT
VEKLLTSESKLKQAALTLEVERSALLQMVEELRRQSARPSTPEPDCTQLEPTGD

Format : Rabbit Whole Serum.

Constitution : 0.05 % NaN₃

Quantity : 100 µL



Specificity : Specific to recombinant protein GX0168. This antibody detects TBC1D1 protein. It also recognizes mouse TBC1D1 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by  ProteinExpress Co., Ltd. and  Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619