



ProteinExpress

ProteinExpress Co., Ltd.
2-1-6, Kazusakamatari,
Kisarazu-shi, Chiba 292-0818, Japan
Tel: +81-438-52-0112, Fax: +81-438-52-0113
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse MACF1 (KIAA1251), Rabbit Polyclonal

Code # : MKA1251 (Lot# xxx)

Gene : Microtubule-actin cross-linking factor 1, isoforms 1/2/3/5 (Actin cross-linking family protein 7) (Macrophin-1) (Trabeculin-alpha) (620 kDa actin-binding protein) (ABP620)

Immunogen : GX2446 (264 amino acids)
ATERYTTLHSHKICIRLGSLSMLLQYQQFQSSADSLQAWVLTCEASVGKLLSDTVASDPG
VLQQQLATTKQLQEELAEHQVPVEKLQKAAHLLDIEGEPALDCRPIQETTDSSSRFQN
LSCSLDERSALLQKAIQSQSVQESMESLLQSIREVEQNLERDQVASLSSGVIQEALANN
MKLKQDIARQKSSLEATHDMVTRFMETADSNASVLOGKLAELSQRFFQQLLQQQEKESN
LKKLLPQAEMFEQLSNKLQQFMEN

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX2446. This antibody detects mouse MACF1 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute

2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>