



# ProteinExpress

ProteinExpress Co., Ltd.  
2-1-6, Kazusakamatari,  
Kisarazu-shi, Chiba 292-0818, Japan  
Tel: +81-438-52-0112, Fax: +81-438-52-0113  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Mouse SDK2 (KIAA1514), Rabbit Polyclonal

**Code # :** MKA1514 (Lot# xxx)

**Gene :** protein sidekick-2 precursor

**Immunogen :** GX2438 (165 amino acids)  
RGQSKKYSKKTDGGNTKSGALGHGEMLSLDESSFPALNNRRLSVKNSFCRKNGLYTR  
SPPRPSPGSLHYSDVDVTKYNDLIPAESSSLTEKPSEISDSQGSSEYEVDTNTQKAHSF  
VNHYISDPTYYNWRRQQKGISRAQAYSYTESSDGEPDHVTVPNS

**Format :** Rabbit Whole Serum

**Preservative :** 0.05 % NaN<sub>3</sub>

**Quantity :** 100 µL

**Specificity :** Specific to recombinant protein GX2438. This antibody detects mouse SDK2 protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at -20°C. Avoid freeze-thaw cycles.

**References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**  
2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>