



ProteinExpress

ProteinExpress Co., Ltd.
2-1-6, Kazusakamatari,
Kisarazu-shi, Chiba 292-0818, Japan
Tel: +81-438-52-0112, Fax: +81-438-52-0113
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti- Mouse METTL14 (KIAA1627) Polyclonal Antibody, Rabbit

Code # : MKA1627PA, (Lot# xxx)

Gene : mouse methyltransferase like 14 (METTL14)

Immunogen : > GX1436 mbg20692 mKIAA1627 167aa(248-748)
EIAETRETCRASYDTSAPNSKRKCLDEGETDEDKVEEYKDELEMQQEEENLPYEEIYKD
SSTFLKGTQSLNPHNDYCQHFVDTGHRPQNFIRDVGLADRFEEYPKLRELIRLKDELIK
SNTPPMYLQADIEAFDIRELTPKFDVILPEPPLEYYRETGITANEK

Format : Rabbit IgG purified with Protein A affinity chromatography.

Preservative : PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity : 100 µg (200 µL)

Specificity : Specific to recombinant protein GX1436. This antibody detects mouse KIAA1627 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>