



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse POLR3H (KIAA1665) Polyclonal Antibody, Rabbit

Code # : MKA1665PA [FIG00883], (Lot# XXX)

Gene : mouse polymerase (RNA) III (DNA directed) polypeptide H (22.9kD)
(mPOLR3H, mKIAA1665)

Immunogen : GX1921 (GST-fusion protein, 107 amino acids)
KGCSPEGVHVS LGFFDILIPPE SLQQPAKFDEAEQVWVWEYETEEGAHDLYMDTGEEIRFRV
VDESFVDTSP TGPSSAEAASSEELPKKEAPYTLVGSISEPGLG

Format : Rabbit IgG purified with Protein A affinity chromatography.

Constitution : PBS containing with 50% glycerol and 0.02% of NaN_3

Quantity : 100 μg (200 μL)

Specificity : Specific to recombinant protein GX1921. This antibody detects mPOLR3H protein. It also recognizes human POLR3H protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C . Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>