



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti- Mouse MTMR12 (KIAA1682) Polyclonal Antibody, Rabbit

Code # : MKA1682 [FSR00265], (Lot# XXX)

Gene : mouse myotubularin related protein 12 (mMTMR12, mKIAA1682)

Immunogen : GX0851 (GST-fusion protein, 247 amino acids)
DTNMGRENLDAQSKPLTLLTVVDWSVQFEPKAQTLLRNPLYVEKPKLDRGQQKGSRFKHQRQ
LSLPLTQSKSSPKRGFFREETDHLIKNLLGKRISKLINSSDDLQDSSREFYDSWHSKPTDYHGLL
LPHIEGPEIKVVAQRYLRWIPEAQILGGGRVATMGKLEMMEEVQSLQEIEARHHRQEAVHAE
APGLLRNSVRLSSLPFALLQRHSAKPVLPSTSGWKALGGEDDLAKREDEFVDLGDV

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX0851. This antibody detects mMTMR12 protein. It also recognizes human MTMR12 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>