



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti- Mouse FAM5B (KIAA1747) Polyclonal Antibody, Rabbit

**Code # :** MKA1747PA [FIG00276], (Lot# XXX)

**Gene :** mouse family with sequence similarity 5, member B (mFAM5B, mKIAA1747)

**Immunogen :** GX0448 (GST-fusion protein, 121 amino acids)  
YEPLEMTDPSKNLGYMKINTLQVFGYSLPFDPAIRDILQLDYPYTGGSQDSALLQLIELRDRV  
NQLSPPGKVRDLDFSCLLRHRLKLANNEVGRIQSSLRAFNSKLPNPVEYETGKLCS

**Format :** Rabbit IgG purified with Protein A affinity chromatography.

**Constitution :** PBS containing with 50% glycerol and 0.02% of NaN<sub>3</sub>

**Quantity :** 100 µg (200 µL)

**Specificity :** Specific to recombinant protein GX0448. This antibody detects mFAM5B protein. It also recognizes human FAM5B protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at -20°C. Avoid freeze-thaw cycles.

**References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**

2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>