



ProteinExpress

ProteinExpress Co., Ltd.
2-1-6, Kazusakamatari,
Kisarazu-shi, Chiba 292-0818, Japan
Tel: +81-438-52-0112, Fax: +81-438-52-0113
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse MYO15B (KIAA1783), Rabbit Polyclonal

Code # : MKA1783 (Lot# CCAI001)

Gene : myosin XVB pseudogene.

Immunogen : GX2418 (246 amino acids)
QILRDTFTESCSTRISQDERHKMKGLLDLEVSLETLDIVEDSIKKRIVVAARDNWANYFS
RIFPVSGESGSDVQLLGVSHRGLRLLKVTQSPSFHLDQLKTLCSYSYAEVLTVQCRGRST
LELSLKNEQLILHTAWARAIAKAMVDLFLSELRKDSGYVIALRSYITDDNSLLSFHRGDLI
RLLPVTTLEPGWQFGSAGGRSGLFPDDVVPAAAPDLSFSLGKRNSWQRKSKLGPAAQEV
KTEEVK

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX2418. This antibody detects mouse MYO15B protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute

2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>