



ProteinExpress

ProteinExpress Co., Ltd.
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
Tel: +81-43-202-5755, Fax: +81-43-202-5756
E-mail: antibody@proteinexpress.co.jp
Web: www.proteinexpress.co.jp

Anti-Mouse EHMT1 (KIAA1876) Polyclonal Antibody, Rabbit

- Code # :** MKA1876AF, (Lot# XXX)
- Gene :** mouse Euchromatic histone-lysine N-methyltransferase 1 (Histone H3-K9 methyltransferase 5, G9a-like protein 1, EHMT1) (mEHMT1, mKIAA1876)
- Immunogen :** GX0523 (GST-fusion protein, 156 amino acids)
CEYVGELISDSEADVREEDSYLFDLDNKDGEVYCIDARFYGNVSRFINHHCEPNLVP
VRVFMHQDLRFPRIAFFSTRLIQAGEQLGFDYGERFWDVKGKLFSCRCGSSKCRH
SSAALAQQRQASAAQEPQENGLPDTSSAAAADPPLILMKCGFAF
- Format :** Affinity Purified Rabbit IgG
- Constitution :** PBS containing with 50% glycerol and 0.02% of NaN_3
- Quantity :** 50 μg (250 μL)
- Specificity :** Specific to recombinant protein GX0523. This antibody detects mEHMT1 protein. Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Other applications have not been tested.
- Storage :** Store at -20°C . Avoid freeze-thaw cycles.
- References :** 1) Okazaki, N. *et al.*: *DNA Res.*, 10(1), 35 (2003).

For research use only, not for diagnostic use



Kazusa DNA Research Institute
2-6-7, Kazusa-kamatari,
Kisarazu, Chiba 292-0818 JAPAN
<http://www.kazusa.or.jp>