



Anti- Mouse RGNEF (KIAA1998) Polyclonal Antibody, Rabbit

Product Code: MKA1998

Gene: mouse Rho-guanine nucleotide exchange factor (mRGNEF, mKIAA1998)

Immunogen: GX1146 (GST-fusion protein, 194 amino acids)

RQEKDVIKRQDVIFELMQTEVHHIQTLLIMSEVFRKGMKEELQLDHSTVDKIFPCLDELLETHRHF
FFSMKERRQESCAGSDRNFVINQIGDILVQQFSEENASKMKRIYGEFCSHHKEAMSLFKELQQN
KKFQNFIKIRNSNLLARRRGIPECILLVTQRITKYPVLVERILQYTKERTEEHRDLCKALGLIK
ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA1998](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL



Specificity: Specific to recombinant protein GX1146. This antibody detects mRGNEF protein. It also recognizes human RGNEF protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki,N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by  ProteinExpress Co., Ltd. and  Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619